Starting Insulin – a patient guide

INSULIN SYRINGES AND PENS

Insulin is injected in the fat just under the skin, using:

- Syringes
- Insulin pens
- Insulin pumps

The most common way to inject insulin is with a syringe.

- A syringe is a hollow plastic tube with a plunger inside and a short skinny needle attached.
- Insulin is injected into the fatty tissue just under the skin. This is called a subcutaneous tissue, or "sub-Q" injection.



Syringes come in different sizes.

- Each line on a 100-unit syringe marks 2 units of insulin.
- Each line on a 50-unit or 30-unit syringe marks 1 unit of insulin.
- Use a syringe large enough to hold the whole dose of insulin.

Place label here.

The Nurse Practitioner Healthcare Foundation (NPHF) acknowledges the Association of Clinicians for the Underserved for their expertise in the initial development of this public education series, Starting Insulin. We also thank Nipro Diagnostics for providing initial funding. Funding for this 2020 update provided by NPHF.



Starting Insulin – a patient guide

INSULIN SYRINGES AND PENS

- Use a 30-unit syringe if you take 30 units of insulin or less.
- Use a 50-unit syringe if you take 50 units of insulin or less.
- Use a syringe that shows 1/2-unit marks if you need 1/2 a unit of insulin.
- Be sure that you can clearly see the markings on your syringe.
- No prescription is needed for insulin syringes.
- If you have poor eyesight or arthritis in your hands, talk to your healthcare provider about using another method, such as an insulin pen.
- Your healthcare provider can show you the different sizes of syringes and help you choose what works best for you.

Needles are described by length and thickness ("gauge").

- The standard needle is 1/2-inch long.
- Needles also come in 5/16-inch and 3/16-inch lengths.
- The 3/16-inch length is often used for children.
- The thinner the needle, the higher its gauge. For example, a 31-gauge needle is thinner than a 28-gauge needle.

Insulin pens look like writing pens, except that there is a thin, short needle at the end.

- Some insulin pens can be refilled, while other pens are thrown away when empty.
- Pre-filled insulin pens come with either one type of insulin or a mixture of two types of insulin.
- Insulin pens with pre-mixes work if they match your prescription.
- You may need one insulin pen for each type of insulin if pre-mix does not match your prescription.

Insulin pumps are used by people who have type 1 diabetes. People with type 2 diabetes rarely use an insulin pump. Insulin pumps give a continuous dose of insulin. Talk to your healthcare provider if you think an insulin pump might be right for you.

